
Operator's Guide

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English:

An issue or revision date for this manual is shown on the front cover.

If you purchased your unit prior to the date listed on the front cover of this manual, the features and operations may not match with those contained in this manual.

ZOLL, M Series, and R Series are trademarks or registered trademarks of ZOLL Medical Corporation in the U.S. and/or other countries. All other trademarks and registered trademarks are the property of their respective owners.

Copyright 2016 ZOLL Medical Corporation. All rights reserved.

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0123

Table of Contents

Degree of Protection Against Electrical Shock.....	1
Symbols Used on the Equipment.....	1
Warranty.....	1
Service.....	1
ZOLL Autoclavable Internal Handle Sets and Internal Paddles	2
Internal Paddle Installation	5
Cleaning and Sterilization	6
Verification of Operation Prior to Use.....	9
Defibrillation Procedure	10
Three Month Checkout Procedure	11
Date Code Identification	12
Ordering Additional Parts	13

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
Degree of Protection Against Electrical Shock

The ZOLL Internal Handles and Paddles are classified as Type CF, Defibrillator Protected equipment.

Symbols Used on the Equipment

Any or all of the following symbols may be used in this manual or on this equipment:

-  Type B patient connection.
-  Type BF patient connection.
-  Type CF patient connection.
-  Defibrillation protected Type BF patient connection
-  Defibrillation protected Type CF patient connection
-  **ATTENTION** Refer to manual for more information
-  Fusible link.
-  Protective (earth) ground terminal
-  Equipotentiality
-  **DANGER** High voltage present
-  Alternating current.
-  **DANGER** Risk of explosion if used in the presence of flammable anesthetics
-  **Conformité Européenne** Complies with the Medical Device Directive 93/42/EEC.

Warranty

(a) ZOLL Medical Corporation warrants to the Customer that from the date of installation or thirty (30) days after the date of shipment from ZOLL Medical Corporation's facility, whichever first occurs, the Internal Handle Sets and Internal Paddles will be free from defects in material and workmanship under normal use and service for the period of 90 days. During such period ZOLL Medical Corporation will, at no charge to the customer, either repair or replace (at ZOLL Medical Corporation's sole option) any part of the equipment found by ZOLL Medical Corporation to be

defective in material or workmanship. If ZOLL Medical Corporation's inspection detects no defects in material or workmanship, ZOLL Medical Corporation's regular service charges shall apply. (b) ZOLL Medical Corporation shall not be responsible for any equipment defect, the failure of the equipment to perform any function, or any other nonconformance of the equipment, caused by or attributable to: (i) any modification of the equipment by the customer, unless such modification is made with the prior written approval of ZOLL Medical Corporation; (ii) the use of the equipment with any associated or complementary equipment, accessory or software not supplied by ZOLL Medical Corporation (iii) any misuse or abuse of the equipment; (iv) exposure of the equipment to conditions beyond the environmental, power or operating constraints specified by ZOLL Medical Corporation; or (v) installation or wiring of the equipment other than in accordance with ZOLL Medical Corporation's instructions. (c) This warranty does not cover items subject to normal wear and burnout during use, including but not limited to lamps, fuses, batteries, patient cables and accessories. (d) The foregoing warranty does not apply to software included as part of the equipment (including software embodied in read-only memory, known as "firmware"). (e) The foregoing warranty constitutes the exclusive remedy of the customer and the exclusive liability of ZOLL Medical Corporation for any breach of any warranty related to the equipment supplied hereunder. THE WARRANTY SET FORTH HEREIN IS EXCLUSIVE AND ZOLL MEDICAL CORPORATION EXPRESSLY DISCLAIMS ALL OTHER WARRANTIES WHETHER WRITTEN, ORAL, IMPLIED, OR STATUTORY, INCLUDING BUT NOT LIMITED TO ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.

For additional information, please call ZOLL Medical Corporation at 1-800-348-9011 (in Massachusetts: 1-978-421-9655). International customers should call the nearest authorized ZOLL Medical Corporation service center.

Service

U.S.A. Customers

Should the ZOLL Internal Handle Sets or Internal Paddles require service, contact ZOLL Medical Corporation at 1-800- 348-9011 (in Massachusetts: 1-978-421-9655).

You will be given a Return Claim number under which to return the ZOLL Internal Handle Sets or Internal Paddles. Send the unit in its original packaging to:

ZOLL Medical Corporation
269 Mill Rd.
Chelmsford, MA 01824-4105
Attn: Technical Service Department

Loaner instruments are available for use while repairs are being completed by ZOLL's Technical Service Department.

For additional information on the product, its preparation for use, cleaning, sanitizing, sterilizing or other questions on infection control procedures for this product, please contact ZOLL Technical Service at 1-800-348-9011 (in U.S.A. and Canada)

International Customers

Should the ZOLL Internal Handle Sets or Internal Paddles require service, they should be returned, in their original container or equivalent, to the nearest authorized ZOLL Medical Corporation service center.

ZOLL® Autoclavable Internal Handle Sets and Internal Paddles

General

The ZOLL Autoclavable Internal Handles are designed for use with the ZOLL **R Series**® and **M Series**® devices and internal defibrillation paddles to defibrillate the heart during open chest surgical procedures.

The Operator's Guide for your ZOLL Defibrillator must be used in conjunction with this guide.

Notes:

- This manual describes the use of the ZOLL Autoclavable Internal Handle Sets only. For instructions related to the use of other ZOLL Internal Handles, see the applicable users guide.
- Cleaning and sterilization protocols for other ZOLL Internal Handle Sets differ significantly and the proper Operator's Guide must be followed.
- ZOLL's Internal Paddles (various sizes) are compatible with all ZOLL Internal Handle Sets.
- The ZOLL Autoclavable Internal Handle Sets are intended for use by or under the direction of a physician.

M Series units only

DO NOT OPERATE AUTOCLAVABLE INTERNAL HANDLES UNLESS YOUR M SERIES DEVICE HAS SOFTWARE REVISION 3.5 OR GREATER.


To determine the revision of your **M Series** device, perform the following steps :

1. Ensure the **M Series** device has been turned OFF for more than 10 seconds.
2. Turn the Selector Switch to **MONITOR**.
3. Observe the system software revision, displayed briefly in the lower right corner of the screen (during power up).
4. If the system software revision is less than 3.5, call ZOLL Technical Service Department at 1-800-348-9011 to request a software upgrade kit.


M Series devices equipped with an advisory ECG analysis feature will only operate in the manual Defibrillation mode when the ZOLL Autoclavable Internal Handles are attached.

ZOLL AUTOCLAVABLE INTERNAL HANDLES AND PADDLES

Autoclavable Internal Handles are available in two styles (shown below). **Style A** devices (P/N 8011-0501-xx) include a button on the apex handle whose activation initiates discharge of the attached defibrillator. **Style B** devices (P/N 8011-0500-xx) include no discharge button. Discharge of the defibrillator is affected by depression of the front panel “SHOCK” button when these devices are used (R Series unit shown below).


Style A : Part # 8011-0501-xx


Style B : Part # 8011-0500-xx

OPERATOR'S GUIDE

All of ZOLL's Internal Handle Sets can be used with interchangeable adult or pediatric Internal Paddles of various sizes. Selection of the appropriate Internal Paddle **diameter** and **shaft length** for therapeutic purposes is the responsibility of the user.

ZOLL Internal Paddle diameters are shown below in approximate scale for reference only.
(See *Ordering Additional Parts* at end of this manual for part numbers and shaft lengths.)

Connection of the Internal Handle Set to the appropriate defibrillator automatically causes the defibrillator to limit its energy output to a maximum of 50 Joules.

Internal Paddle Installation

The Internal Handle Sets and Internal Paddles can be sterilized by either steam sterilization in a pre-vacuum autoclave according to the guidelines that we provide in this manual or by the STERRAD® Sterilization System.

ZOLL's Internal Handle Sets are sealed to withstand up to 100 steam sterilization cycles or 100 STERRAD® sterilization cycles, as specified by the following table:


STERRAD® Sterilization Systems	Cycle Selection
STERRAD® 50	N/A
STERRAD® 100S and 200	Short
STERRAD® NX®	Standard
STERRAD® 100NX™	Standard

For details on cycle selection definitions, refer to the STERRAD® System User's Guide.


Do not disassemble the Internal Handle Sets. Any attempt to disassemble the Internal Handle Sets will void any applicable warranties.

To attach the paddle to the handle set:

1. Carefully thread each Internal Paddle into the Internal Handle.
2. Tighten paddles to the maximum number of turns allowed by the Handle Assembly, until the Paddle's flange and Handle top meet together snugly.

3. (*with Discharge button only*) Unscrew paddle up to one full turn to set the Spoon Face and the Discharge button into an appropriate orientation for use.

CAUTION

- Do **not** over tighten the Internal Paddles.
- Remove the Internal Paddles from service after 100 uses. Normal wear of the Paddle's threads may degrade the electrical contact between the Paddle and the Hand Assembly.

Cleaning and Sterilization

Notes:

- The following process for cleaning and sterilization of the Autoclavable Internal Handle Sets and Internal Paddles has been validated to be effective in the disinfection of these products. Product users are responsible for qualifying any deviation from this recommended method of processing.
- All cleaning agents should be prepared and used according to the manufacturer's instructions and product labeling.

Do **not** allow the Internal Handle Sets or Internal Paddles to dry after use and before cleaning. Items contaminated with blood and/or other protein materials cannot be effectively cleaned if allowed to dry.

Internal Handle Sets (including the connector) and Internal Paddles may be immersed, if necessary during cleaning.

Internal Paddles must be disassembled from the Internal Handles prior to cleaning.

The Internal Handle Sets and Internal Paddles must be hand washed. Do **not** subject them to ultrasonic or mechanical washing machines.

To ensure the integrity of the Internal Handle Sets and Internal Paddles, perform a functional test of the assembled device prior to each use (see "Verification of Operation Prior to Use").

Internal Handle Sets and Paddles

CAUTION

- The ZOLL Internal Handle Sets and Internal Paddles are sold and shipped in a non-sterile condition. They should be cleaned and sterilized following the procedures outlined below prior to their first use and after each reuse.
- Do not drop, bump or knock the Internal Paddles. Damage to the insulating coating of the Internal Paddle may occur.
- Inspect each Internal Paddle after cleaning and prior to each use for:
 - Nicks or burrs that may injure patient tissue;
 - Scratches, pits, or gouges in the Internal Paddle surface;
 - Loose, flaking or damaged insulating coating.If any of these conditions are observed, remove the Internal Paddle from use.
- Internal Paddles must be removed from handles prior to cleaning and sterilization.
- ZOLL Internal Handle Sets are constructed from high quality materials. The severe conditions of sterilization, however, will limit their useful life. Thus, the useful life of the Internal Handle Set is limited primarily by the frequency of sterilization rather than the age of the Internal Handle Set.

- Inspect the Internal Handle Set frequently for signs of deterioration such as cracks, crazing, damaged cables, and damaged switch covers. Replace if deterioration is noted.
- Do not clean the Internal Handle Sets in mechanical washers or ultrasonic cleaners. Hand wash only.
- Do not expose the Internal Handle Set to any product containing organic solvents such as acetone, ketones, chlorinated hydrocarbons or aromatic hydrocarbons. Exposure to these solvents may degrade the handle materials.
- Keep the connector as free from contamination as possible. The connector may be immersed and cleaned with the rest of the assembly but is difficult to clean properly because of its shape and function.
- The Connector cap must remain in place at all times to protect the contact sockets, except when in use and during cleaning.
- Repeated sterilization cycles may cause the cable insulation to blister or bubble. This is considered normal and will not affect the function of the handle set. Replace the handle set if the cable insulation is cracked or cut.

1. Disassemble

Disassemble Internal Paddles from the Internal Handle Set prior to cleaning by turning paddles counterclockwise.

2. Wash

Internal Handle Sets, Paddles, cable connector and connector cap must be thoroughly cleaned prior to each sterilization.

Prepare enzymatic detergent as recommended by the manufacturer's product labeling (e.g. solution temperature, concentration, etc.). Submerge the Handle Set or Paddles in enzymatic detergent and scrub with a soft bristled brush. Soak for 5 minutes.

Note: It is recommended that a brush with a ½ inch (1.2 cm) diameter be used for cleaning the internal cavities.

Following the soak period, scrub the Handle Set, Paddles, cable connector and connector cap again with a soft bristled brush. Visually inspect the Handle Set, Paddles, cable connector and connector cap for cleanliness. If necessary, repeat the process above. The Handle Set, Paddles, cable connector and connector cap must be free of any visible contamination prior to rinsing.

3. Rinse

Rinse the Handle Set, Paddles, cable connector and connector cap for at least 30 seconds under running water. When rinsing the Internal Handle Set, be sure to rinse the cable connector, connector cap and the handle end thoroughly.

If the manufacturer of the cleaning agent recommends a rinse protocol, this should be followed.

4. Air Dry

Allow the Handle Set, Paddles, cable connector and connector cap to air dry prior to wrapping for sterilization. Air drying may be performed in a drying oven whose temperature does not exceed 120° C (250° F).

5. Inspect

Inspect the clean Handle Set, Paddles, cable connector and connector cap for any residual contaminants, particularly in the area where the Internal Paddles thread into the handles.

Inspect Internal Handle Sets for signs of deterioration such as cracks, crazing, damaged cables, connector pins or switch covers.

Repeated sterilization cycles may cause the cable insulation to blister or bubble. This is considered normal and will not affect the function of the handle set. Replace the handle set if the cable insulation is cracked or cut.

Inspect Internal Paddles for loose, flaking or damaged insulating coating, or other mechanical damage such as: scratches, pits, gouges, nicks or burrs that may injure patient tissue.

If any of these conditions are observed, remove the Internal Handle Set or Paddles from use.

6. Secure the Connector Cap

Attach the connector cap, located on the cable of the Handle Set, securely to the cable connector (Refer to **Figure 1**).

The Connector cap must remain in place at all times to protect the contact sockets, except when in use and during cleaning.


Figure 1

7a. Sterilize – Steam Sterilization Method

Wrap: Before wrapping the Internal Handle Set, be sure the cable is coiled with a diameter larger than 6 inches (15 cm). Never wrap the cables around the handles. This may hinder proper sterilization.

Wrap the Internal Handle Set or Paddles in sterile wrapping sheets per your institution's wrapping procedures.

Sterilize: Sterilize the wrapped Internal Handle Set or Paddles using a pre-vacuum steam autoclave. The autoclave cycle should have the following characteristics:

1. Temperature: 132.2° - 135° C (270° - 275° F).
2. Exposure time: 4-5 minutes.
3. Drying time: 10-12 minutes.

The Internal Handle Sets and Internal Paddles can withstand up to 100 steam sterilization cycles.

7b. Sterilize – STERRAD® System Hydrogen Peroxide Gas Plasma Sterilization Method

Wrap: Before wrapping the Internal Handle Set, be sure the cable is coiled with a diameter larger than 6 inches (15 cm). Never wrap the cables around the handles. This may hinder proper sterilization.

Wrap the Internal Handle Set or Paddles in sterile wrapping sheets per your institution's wrapping procedures.

The ZOLL Internal Handle Sets are compatible with the STERRAD® NX, 100S, 100NX, 50, and 200 sterilization equipment systems. The Internal Handle Sets can withstand up to 100 STERRAD® cycles.

Verification of Operation Prior to Use

WARNINGS

- Do not use the ZOLL Autoclavable Internal Paddles in the presence of flammable agents, oxygen rich atmospheres, or flammable anesthetics. Using the paddles in the presence of such agents may cause an explosion.
 - Users of the ZOLL Internal Paddles should assess the defibrillator's performance in their typical environment of use for the possibility of interference from high power radio or electro-surgical units. This interference may be observed as shifts in monitor baseline, trace compression, or transient spikes on the defibrillator display.
 - Verify that no one is in contact with monitoring cables, leads, bed rails, or any other potential current pathway prior to defibrillator discharge.
 - All persons nearby must be warned to “**STAND CLEAR**” prior to defibrillator discharge.
 - When performing pre-use checks, keep hands away from the Internal Paddles while pressing the **DISCHARGE** button.
 - When verifying high voltage wiring, paddles must be held **firmly** together so that they are not damaged.
 - Inspect each Internal Paddle prior to use for:
 - Nicks or burrs that may injure patient tissue;
 - Scratches, pits, or gouges in the Internal Paddle surface;
 - Loose, flaking or damaged insulating coating.
- If any of these conditions are observed, remove the Internal Paddle from use.

- 1.) Connect the Internal Handle Set with attached Internal Paddles to the defibrillator.
- 2.) Verify that the **DISCHARGE** button does not stick. Before charging the defibrillator, press the **DISCHARGE** button and verify that there is an audible click and that the button springs back upon release.
- 3.) Charge the defibrillator to 2 Joules. Wait for the READY tone.
- 4.) During the READY tone, hold both handles and Internal Paddles out and away from any person or object, press and hold the **DISCHARGE** button located on one of the handles to simulate delivery of energy to the patient.
- 5.) Verify that the defibrillator doesn't discharge and displays the error message, “POOR PAD CONTACT”. Occurrence of this message verifies that the “**SHOCK**” button located on the apex handle is operating correctly.
- 6.) To verify the integrity of high voltage wiring contained in the handle set. Charge the unit to 30 Joules then place the paddles firmly together and press the “**SHOCK**” button. The device will discharge and display the messages “TEST OK” and “XX J DELIVERED” (xx = 28-30 Joules).

WARNING: When performing pre-use checks, keep hands away from the Internal Paddles while pressing the DISCHARGE button.

Defibrillation Procedure


WARNINGS

Before proceeding, CAREFULLY read the following:

- Defibrillation should be attempted only by appropriately trained and skilled personnel who are familiar with equipment operation.
- Do not touch the bed, patient, or any equipment connected to the patient during defibrillation.
- All persons in attendance of the patient must be warned to “**STAND CLEAR**” prior to defibrillator discharge.
- Do not touch Internal Paddles together during discharge into a patient.
- Do not detach Internal paddles from the handles assembly in the vicinity of open surgical wounds. Loose particles may be produced if these parts have been cross threaded.

Review your ZOLL defibrillator's Operator's Guide's to become familiar with the defibrillator's operation. Note the following special instructions for use of internal paddles.

- 1) Connect the Internal Handle Set with attached Internal Paddles to the defibrillator.
- 2) Press the up (▲) and down (▼) arrows of the **Energy Select** button located on the front panel of the **M Series** device to select the desired energy.
- 3) Press the **CHARGE** button on the defibrillator front panel to charge the unit to the selected energy level. The maximum selected level is limited to 50 Joules.
- 4) Hold Internal Paddles as shown in following figure. Position the Internal Paddles appropriately on the patient's heart.

- 5) **Style A Handles:** Press and hold the **DISCHARGE** button on the apex handle to deliver the selected energy to the patient.
- 6) **Style B Handles:** Handles without a discharge button on the apex handle, press and hold the **SHOCK** button on the front panel to deliver the selected energy to the patient.

Three Month Checkout Procedure

To ensure quality operation of the Internal Handles perform this Checkout Procedure at least every 3 months.

Continuity


Remove the Internal Handles from the defibrillator for this test.

Using an electrical continuity tester, such as a volt/ohm meter (VOM) or a digital multimeter (DMM), verify the resistance between the test points in the table on the next page.

Be sure that good electrical contact is made between the test probes and the test points.


M Series Internal Handles (No Button on Handle)

Pin 5 to Sternum Paddle	< 1 ohm
Pin 1 to Apex Paddle	< 1 ohm
Pin 2 to Pin 3	Between 1.42 and 1.52 K ohms

M Series Internal Handles (with Button on Handle)

Pin 5 to Sternum Paddle	< 1 ohm
Pin 1 to Apex Paddle	< 1 ohm
Pin 2 to Pin 4 with shock button pressed	Between 2.72 and 2.88 K ohms
Pin 3 to Pin 2	Between 882 and 936 ohms

Leakage Current


Connect the Internal Handles and Paddles to the defibrillator. Perform a standard electrical safety leakage test. The system leakage current should not exceed 50 μA @ 132 VAC or 100 μA @ 264 VAC.

Date Code Identification

To determine the date of manufacture for the Internal Paddles, the date codes can be found on the front and back of the base of the paddle.

Front: The two numbers indicate the year of manufacture (Example: 15=2015, 16=2016, etc.).

Back: The first character is a letter from A-L that indicates the month of manufacture (Example: A=January, B=February and so on to L for December.) The second character is a number from 1-5 that indicates the week of the month (Example: 1=first week, 2=second week and so on to 5 for fifth week.)

Ordering Additional Parts

Reorder numbers for the parts most frequently ordered are listed below:

Part Number	Item Description	Adult	Pediatric
8011-0500-xx	Autoclavable Internal Handles with No Button		
8011-0501-xx	Autoclavable Internal Handles with Button		
8000-0277-01	Overmolded Internal Paddle, 2.7" (6.8 cm) dia. X 5.6" (14.3 cm) Short pair	●	
8000-0277-02	Overmolded Internal Paddle, 2.7" (6.8 cm) dia. X 8.1" (20.6 cm) Long pair	●	
8000-0277-03	Overmolded Internal Paddle, 2.0" (5.1 cm) dia. X 4.9" (12.5 cm) Short pair		●
8000-0277-04	Overmolded Internal Paddle, 2.0" (5.1 cm) dia. X 7.4" (18.7 cm) Long pair		●
8000-0277-05	Overmolded Internal Paddle, 1.6" (4.0 cm) dia. X 4.5" (11.4 cm) Short pair		●
8000-0277-06	Overmolded Internal Paddle, 1.6" (4.0 cm) dia. X 7.0" (17.7 cm) Long pair		●
8000-0277-07	Overmolded Internal Paddle, 1.0" (2.5 cm) dia. X 3.9" (10.0 cm) Short pair		●
8000-0277-08	Overmolded Internal Paddle, 1.0" (2.5 cm) dia. X 6.4" (16.3 cm) Long pair		●
8000-0277-09	Overmolded Internal Paddle, 3.0" (7.6 cm) dia. X 5.9" (15.0 cm) Short pair	●	
8000-0277-10	Overmolded Internal Paddle, 3.0" (7.6 cm) dia. X 8.4" (21.3 cm) Long pair	●	

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