



VS-600 supports fast, accurate and comfortable NIBP measurements. It displays reliable temperature readings in seconds and supports oxygen saturation measurements with your choice of three state-of-the-art SpO<sub>2</sub> technologies. With all these features VS-600 will greatly reduce the time needed for vital signs checking, allowing you to focus even more on your patients.

**mindray**

**VS-600**  
Vital Signs Monitor

Clear, Simple and Portable



[www.mindray.com](http://www.mindray.com)

P/N:ENG-VS-600-210285X6P-20190529  
©2019 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.

 **procamed**  
cardiac medical systems

Procamed AG  
Schulstrasse 5  
CH-8355 Aadorf

Tel+41 52 368 60 00  
[www.procamed.ch](http://www.procamed.ch)  
[info@procamed.ch](mailto:info@procamed.ch)

**mindray**  
healthcare within reach



5.7" large and bright display

Large and easy-to-read numerics


Intuitive user interface design

Designed with an upward inclination view angle


VS-600 is an ideal portable vital signs monitor, satisfying your basic clinical needs for periodic spot-check in the low-acuity areas of your hospital. VS-600 is suitable for adult, pediatric and neonatal patients and features a clear display, simple operations and superior portability.

## Portable





 Supper Li-ion battery for continuous use up to 22 hours




 Compact and light-weight with only 1.7kg (including battery)



 Stable trolley with quick release button



 Robust and comfortable handle and rugged housing



## Functions

**NIBP:**  
Fast and convenient manual mode for spot-check

**Smart Temp:**  
Improved TEMP measurement, first reading within 12s (oral, axillary, and rectal)

**SpO<sub>2</sub>\*:**

- Choices of Mindray, Masimo SET®, Nellcor™ OxiMax™ technologies.
- PI (Perfusion Index) from Mindray and Masimo guides the positioning of SpO<sub>2</sub> measurement

\* Mindray provides 3 options of SpO<sub>2</sub> measurement, Mindray, Masimo and Nellcor. For further information about the availability of Masimo and Nellcor SpO<sub>2</sub>, please contact with your local sales representatives.

## Simple

Easy access to all functions through a four-button keypad

SpO<sub>2</sub> and NIBP accessories compatible with Mindray uMEC, ePM and BeneVision N Series Monitors