

BeneFusion eDS Basic

Infusion Supervision System

Data Sheet



Size		Air pressure	57.0-107.4kPa for operating, 16.0-107.4kPa for storage
1 controller and 4 slots	≤ 270 x 173 x 395 mm	Alarms	
1 controller and 2 slots	≤ 270 x 173 x 245 mm	Type	Audible and visual alarm
Extension 4 slots	≤ 270 x 173 x 335 mm	2 Levels	High and Low, consistent with the highest priority pump alarm.
Weight		Other	
1 controller and 4 slots	≤ 3.4 kg	Mounting type	Horizontal, vertical
1 controller and 2 slots	≤ 2.5 kg	Pole clamp	Horizontal: Diameter range 10 × 25 mm Vertical: Diameter range 15-40mm One pole clamp per shelf module
Extension 4 slots	≤ 3.0 kg	Maintenance Classification	Support 2-16 pumps tool-free assembly Type CF, class I, IP33
Power supply			
AC power	100-240V~, 50/60Hz, 8-3.4A		
Environment			
Temperature	5-40°C for operating, -30-70°C for storage		
Humidity	15-95% for operating, 10-95% for storage		

Mindray Building, Keji 12th Road South,
High-tech Industrial Park, Nanshan, Shenzhen 518057, P.R. China
Tel: +86 755 8188 8998 Fax: +86 755 26582680
E-mail: intl-market@mindray.com www.mindray.com

mindray | Healthcare with reach are registered trademarks or trademarks owned by Shenzhen Mindray Bio-medical Electronics Co., LTD
© 2020 Shenzhen Mindray Bio-medical Electronics Co., Ltd. All rights reserved. Specifications subject to changes without prior notice.
P/N: ENG-Benefusion eDS Basic Datasheet-210285x2P-20201125

mindray